

Appendix F - Blood Analyzers

To add blood analyzers to your organization:

1. Login as an administrators
2. Click on your organization name
3. In the "Blood analyzers" tab
 1. Click add analyzer
 2. Select one of the available option



4. Enter a start date
5. End date is optional and when empty it means that the analyzers is still being used
6. Save

To remove an analyzer from your lab:

1. Login as an administrators
2. Click on your organization name
3. In the "Blood analyzers" tab
 1. Enter an end date to the relevant analyzer
 2. Save: depending on the end date, the analyzer will be moved to the "Blood Analyzers History" session.

Business rules:

- the same blood analyzer cannot be added multiple times for the same period