

Appendix F - Blood Analyzers

To add blood analyzers to your organization:

1. Login as an administrators
2. Click on your organization name
3. In the "Blood analyzers" tab
 1. Click add analyzer
 2. Select one of the available option



The screenshot displays the 'Blood Analyzers' management interface. At the top, there are navigation tabs: 'Organization roles', 'User Rights', 'Lab Info', 'Business rules', and 'Blood Analyzers'. A red arrow points to a yellow 'Add Blood Analyzer' button in the top right corner. Below the navigation, the 'Active Blood Analyzers' section contains a form with the following fields: 'Analyzer*' (text input with 'Sysmex-XT'), 'Start Date*' (calendar icon, '01-Oct-2007'), and 'End Date:' (calendar icon, 'dd-MM-yyyy'). A 'Blood Analyzers History' section is visible below the active analyzers.

4. Enter a start date
5. End date is optional and when empty it means that the analyzers is still being used
6. Save

To remove an analyzer from your lab:

1. Login as an administrators
2. Click on your organization name
3. In the "Blood analyzers" tab
 1. Enter an end date to the relevant analyzer
 2. Save: depending on the end date, the analyzer will be moved to the "Blood Analyzers History" session.

Business rules:

- the same blood analyzer cannot be added multiple times for the same period