

# Notifications (Target: ADOs, Labs, APMUs)

When the uploaded temp logger file contains an error, the Testing Authority (TA) and the APMU of the related sample will be notified. The following errors can generate a notification:

- If there are any temperature value less than 0°C
- All temperature values are  $\geq 0$ , but any value is greater than 50°C.
- If at any point, there isn't one line per minute.
- All values are within 0-50°C, but the first recorded temperature time is after or equal to the sample collection date and time.
- All values are within 0-50°C, but the data logger stop date and time (last recorded date time temperature) precedes the time of reception".

The notification will be sent with the following details:

- **Sender:** the user who generates the activity that triggered the notification
- **Recipients:** TA and APMU users who have their account configured to receive the notification type "New Data Logger Uploaded".
- **Reference Link:** Link to the DCF
- **Content:** The description of the BSS status or uploaded temp logger file

To receive the notifications, follow the steps below:

1. Login with the administrator account;
2. Select "User account Management";
3. Search for the relevant account and go in the "user account info";
4. Add the notification "New Data Logger Uploaded" and save.